
 Ningbo Xinwell Medical Technology Co., Ltd.	Document No.	RD/S-109-0002
	Issued Date	2021.4.15
	Version	A/0

Instructions for Use


Disposable Core Biopsy Instruments

Change History				
Ver.	Change Descriptions	Author	Date	Change No.
1.0	manuscript	Nongke Pan	2020.6.12	/
A/0	A. Add Indications for Use B. Add Intended users C. Add Intended patient group D. Added to the packaging label: Single sterile barrier system with protective packaging outside, Model number, Unique device identifier ,Humidity limitation ,Does not contain natural rubber latex, Non-pyrogenic, Medical device, Quantity E. Reduced in packaging labeling: Gauge size x Needle length(Maatmaat x Naald lengte) ,Fragile, handle with care, Breekbaar - voorzichtig behandelen, Do not roll Niet rollen, Temperaturelimit, Temperatuurgrenzen, Language: English, Taal: Engelse taal, Language: Nederlands, Taal: Nederlands	Nongke Pan	2021.4.15	
Approval Column For Document				
Prepared By		Date	2021.4.15	
Checked By		Date	2021.4.15	
Approved By		Date	2021.4.15	

Instructions for Use

Device Description

Xinwell Core Biopsy Instrument is a single use core biopsy device. It is available in several needle gauge sizes and lengths. The side and rear actuator buttons are color coded according to the various gauge sizes, e.g., Yellow=20G, Pink=18G, Purple=16G and Green=14.

Model	Gauge Size and Needle Length	Length of Sample Notch	Penetration Depth	Button color
BN-OCR-1/1410	14G(2.1mm) ×10cm	19mm	22mm	Green
BN-OCR-1/1416	14G(2.1mm) ×16cm	19mm	22mm	
BN-OCR-1/1610	16G(1.6mm) ×10cm	19mm	22mm	Purple
BN-OCR-1/1616	16G(1.6mm) ×16cm	19mm	22mm	
BN-OCR-1/1810	18G(1.2mm) ×10cm	18mm	22mm	Pink
BN-OCR-1/1816	18G(1.2mm) ×16cm	18mm	22mm	
BN-OCR-1/1820	18G(1.2mm) ×20cm	18mm	22mm	
BN-OCR-1/1825	18G(1.2mm) ×25cm	18mm	22mm	
BN-OCR-1/2010	20G(0.9mm) ×10cm	18mm	22mm	Yellow
BN-OCR-1/2016	20G(0.9mm) ×16cm	18mm	22mm	
BN-OCR-1/2020	20G(0.9mm) ×20cm	18mm	22mm	

Intended use

The disposable core biopsy instruments is intended for use in obtaining biopsies from soft tissues such as liver, kidney, prostate, spleen, lymph nodes and various soft tissue tumors. It is not intended for use in bone.

Indications for Use

Used for obtaining biopsies form soft tissues such as liver, kidney, prostate, spleen, lymph nodes and various soft tissue tumors.

Intended users

The trained medical professionals.

Intended patient group

Unlimited age

Performance characteristics

1. Appearance and cleanliness

- 1) There should be no scars, burrs or foreign material attached to the needle cannula, and the external surface should be smooth without defects.
 - 2) The connection between the needle base and the needle cannula should be straight without obvious skew.
 - 3) Plastic components should be no obvious impurities, burrs, pores and other defects; Handle Logo: unlock, lock symbols should be clear.
2. Needle tip
- Needle stylet and needle cannula edge should be sharp, no visible hook, flat head, burr and other defects.
3. Flowing
- The needle stylet should slide smoothly in the needle cannula without sticking.
4. Handle
- 1) Trigger device
When the needle cannula and stylet are pulled to the limited position, trigger the "pressing button" and "pushing button" respectively for 10 times, the product shall be smooth without jamming.
 - 2) Safety switch
It shall not be activated when the "pushing button" or "pressing button" is triggered when the safety lock is locked.
- 5 Chemical performance
- 1) Ethylene oxide residue: shall be $\leq 4\text{mg/device}$.
 - 2) Corrosion resistance: The part of needle tube and needle core being soaked may not have corrosion trace. The scale mark shall not drop, be missing or illegible.
- 6 Biological performance
- 1) Sterility: After sterilized by ethylene oxide, the product shall be sterile.
 - 2) Bacterial endotoxin: The bacterial endotoxin limit shall be less than 20EU/piece.
 - 3) Pyrogen: There shall be no pyrogen reaction.

Contraindications

1. Medical judgment should be exercised in deciding whether a biopsy should be performed in a patient receiving anticoagulant therapy or in a patient with a bleeding disorder;
2. Infection in the skin at the puncture site.

Warnings

1. Do not use if the package is damaged and the validity period is exceeded!
2. Products are sterilized by ethylene oxide and are non-pyrogenic. For single use only, do not reuse. Reuse of this medical device may present a risk of cross-patient infection because medical devices, especially those with small lumens, interfaces, and/or crevices between components, are difficult or impossible to clean after contact with body fluids or tissues with potential pyrogenic or microbial contamination. Residues of biological agents and contamination with pyrogens or microorganisms leading to infectious complications.
3. Do not resterilize the needle. Sterility of product is not guaranteed after repeated sterilization because the degree of potential pyrogenicity or microbial contamination that may result cannot be determined. Cleaning, reprocessing, and/or resterilization of

medical devices increase the probability of device failure because components may be adversely affected by thermal and/or mechanical changes.

4. The care of patients after biopsy depends on the biopsy method used and the physiological conditions of each patient. Patients should be observed for important physiological parameters and other precautions should be observed to avoid and/or treat potential complications associated with biopsy.
5. Cancerous tissue was found by pathological analysis after biopsy. A "negative" biopsy result in a suspected radiotherapy result does not exclude the possibility of carcinogenesis.
6. The used products shall be isolated and managed as wastes in a centralized manner, and shall not be stored at will. They shall be destroyed and disposed of in accordance with the relevant national laws and regulations, and shall not be arbitrarily discarded.

Precautions

1. This instrument should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied. According to the actual lesion characteristics, the physician selects the imaging guidance method, patient position, sampling method, needle selection, biopsy needle approach, needle insertion method, etc.
2. The introduction of the needle into the body should be carried out under imaging control (ultrasound, X-Ray, CT, etc.).
3. Never test the product by firing into the air. Damage may occur to the needle/cannula tip and could result in patient and/or user injury.
4. Unusual force applied to the stylet or unusual resistance against the stylet while extended out of the supportive cannula may cause the stylet to bend at the specimen notch. A bent specimen notch may interfere with the needle function.
5. Needles should be inspected after each sample collection for needle damage, bent needle cannula, or other defects during collection of multiple samples from the same site. If a defect is found, do not use the associated needle. (No more than 5 sample collections)
6. Please read the instruction carefully before use; Refer to the device catalogue number to determine the device specifications.
7. To avoid infection, asepsis must be maintained throughout the operation.
8. The instrument should be kept away from exposure to high temperature.
- 9 There were no adverse events due to product performance issues.

Potential complications

Potential complications associated with core biopsy procedures are site specific and include, but are not limited to: hematoma; hemorrhage; infection; adjacent tissue injury; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and air embolism. Air embolism is a rare but serious potential complication of lung biopsy procedures. Rapid deterioration of neurological status and/or cardiac arrhythmia may be indicative of air embolism. Prompt diagnosis and treatment must be considered if the patient exhibits signs or symptoms of air embolism.

Directions for Use

The following application methods are the application methods of this product. Other techniques (use of imaging equipment, local anesthesia, small incision, etc.) not described in the biopsy process should be performed based on the experience of all doctors.

Note: The imaging equipment includes B ultrasound, X-ray and CT. Since the needle cannula and needle stylet belong to metal parts, they can be fully developed under the working state of the above imaging equipment. In the operation under imaging equipment, the patient shall cooperate with the doctor's prompt to adjust respiration. In order to improve the accuracy of puncture, it is recommended to use the corresponding disposable puncture stent.

Equipment used in coordination:


- Applicable imaging equipment
- Surgical gloves and surgical drapes
- Local anesthesia
- Scalpel
- Sample collection container

Other necessary equipment


[Biopsy before]

Step 1: Check the biopsy instrument blister package ensuring that there is no damage, remove the biopsy instrument from the blister package, remove the protective sleeve.


Step 2: Set the safety switch to the "  " locking indication side.

Step 3: The needle is in the state shown in the Figure: Check the biopsy instruments whether there is a needle damage, curvature of the needle, or other defects that may affect normal function. Do not use the biopsy needle if it is damaged or bent.

Step 4: Press the button of needle cannula to the bottom, and then press the button of needle stylet to the bottom, and the biopsy needle is in standby state.

Recommendation: For ease of insertion, a scalpel may be used to open the skin at the entry site.

[Biopsy method]

Sterile technique must be used during biopsy.


Step 1: Prepare as needed. Adequate anesthesia should be administered prior to skin incision.

Step 2: Check to confirm that the device is available.


Caution: Do not place your finger on the needle stylet button or the needle cannula button while the device is in standby mode. The functionality of the device is affected when the stylet key and the needle cannula key are blocked.

Step 3: Insert the needle into the lesion edge under the guidance of the imaging device.

Step 4: Push the safety switch to the "🔓" unlocking indication side.


Step 5: While maintaining device position and orientation, press the firing button or push to biopsy tissue.

[Remove specimen]


Step 1: After the biopsy is completed, pull out the biopsy instruments from the patient's body, press the needle cannula button to the bottom again, withdraw the needle cannula, so that the needle tip is exposed to the tissue.

Step 2: Remove the specimen from the needle core sampling groove.

Step 3: To prepare for the next biopsy, press the needle stylet button again to the bottom and

shift the safety switch to the "  " blocking indication side.

Step 4: If the same organ needs to be biopsied again, start again from [Biopsy before] step3 to [Remove specimen] step 4.

Storage and handling condition

During transportation, the instruments shall be protected from heavy weight, direct sunlight and rain, or according to the provisions of the order contract.

The device shall be stored in a clean room where it is protected from light, rodents, fire and insects, well ventilated and without corrosive gas.


Transportation/storage conditions: Relative humidity: $\leq 80\%$

Shelf Life

Sterilized using Ethylene Oxide. Shelf life is 2 years. Manufacture date and expiration date on the product label.

Packaging Marks

Symbol	Title or Description of Symbol	Symbol	Title or Description of Symbol
	Manufacturer		Authorized representative in the European Community
	Country of manufacture: CHINA		Use-by date
	Batch code		Catalogue number

	Sterilized using ethylene oxide		Single sterile barrier system with protective packaging outside
	Do not re-sterilize		Do not use if package is damaged and consult instructions for use
	Model number		Keep away from sunlight
	Keep dry		Unique device identifier
	Humidity limitation		Do not re-use
	Consult instructions for use or consult electronic instructions for use		Caution
	Does not contain natural rubber latex		Non-pyrogenic
	Medical device		Quantity
	The product has passed CE certification		This way up

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Lotus NL B.V.