



Breathing Freely
Breathing Free

Closed suction catheters

A Xinwell Member Company

HUAKUN

Features

1 Safe and high-efficiency

Good air-tightness makes double protection for patients and medical staff, reducing the risk of nosocomial infection; Keeping ventilation during suction to maintain the patient's blood oxygen; the end of the suction catheter is soft, and the mouth is smooth, which does not damage the airway mucosa; repeated suction, decreasing medical waste and reducing medical costs.

2 Easy operation

Ergonomic suction switch design, with locking device; clear and transparent protective cover for easy observation; double screw joint design, easy to move and preventing damage. Equipped with extendable tubes and separators for easy connection and separation.

3 Complete models

For adults and children; 24 hours and 72 hours.

4 Easy to distinguish

The product conforms to the international color code distinguishing model; equipped with a color code card, it is convenient for nursing and reminding replacement information.



Specification

Conventional configuration: ①Extendable tube ②Separator ③Color code card (Date) ④Endotracheal tube conversion (For Children) ⑤Medicine injection kit (For Children)

	NO.	Size (Fr)	Color code	Length (cm)	Type
	H106400 (Children)	6	Light Green	40	H1 (normal type) 24hours
	H108400 (Children)	8	Light Blue	40	H1 (normal type) 24hours
	H110400	10	Black	40	H1 (normal type) 24hours
	H112400	12	White	40	H1 (normal type) 24hours
	H114400	14	Green	40	H1 (normal type) 24hours
	H116400	16	Orange	40	H1 (normal type) 24hours
	H118400	18	Red	40	H1 (normal type) 24hours
	H110600	10	Black	60	H1 (normal type) 24hours
	H112600	12	White	60	H1 (normal type) 24hours
	H114600	14	Green	60	H1 (normal type) 24hours
	H116600	16	Orange	60	H1 (normal type) 24hours
H118600	18	Red	60	H1 (normal type) 24hours	
	H206400 (Children)	6	Light Green	40	H2 (valve type) 72hours
	H208400 (Children)	8	Light Blue	40	H2 (valve type) 72hours
	H210400	10	Black	40	H2 (valve type) 72hours
	H212400	12	White	40	H2 (valve type) 72hours
	H214400	14	Green	40	H2 (valve type) 72hours
	H216400	16	Orange	40	H2 (valve type) 72hours
	H218400	18	Red	40	H2 (valve type) 72hours
	H210600	10	Black	60	H2 (valve type) 72hours
	H212600	12	White	60	H2 (valve type) 72hours
	H214600	14	Green	60	H2 (valve type) 72hours
	H216600	16	Orange	60	H2 (valve type) 72hours
H218600	18	Red	60	H2 (valve type) 72hours	

Ningbo Huakun Medical Equipment Co., Ltd.

Address: No.208, Binjiang RD, High Tech District, Cixi, Zhejiang
 Tel: 86-574-63029812
 Fax: 86-574-63073411
 E-mail: huakun@xinhaigroup.com
 Web: www.chinahuakun.cn

A Xinwell Member Company

